

Torque Setting Reference Chart

Abutment Screw Torque and Driver

Manufacturer	Torque N-cm						Driver	
	15	20	24	25	30	35	Titanium Screw	Gold Screw
Astra 3.0	•						0.050" (1.25 mm) hex	
Astra 3.5 - 4.0		•					0.050" (1.25 mm) hex	
Astra 4.5 - 5.0				•			0.050" (1.25 mm) hex	
BioHorizons					•		0.050" (1.25 mm) hex	
Biomet 3i™ External Hex						•	0.048" hex	Square
Biomet 3i™ Certain® *		•					0.048" hex	
Brånemark (Nobel Biocare)						•	Unigrip	Square
Dentsply Ankylos® C/X	•						1 mm hex	
HiOssen HG Mini (3.5)		•					1.2 mm hex	
HiOssen HG Standard (4, 4.5, 5)					•		1.2 mm hex	
Implant Direct™						•	0.050" (1.25 mm) hex	
MIS® Biocom & Seven						•	0.050" (1.25 mm) hex	
NobelActive™ & Replace™						•	Unigrip	
NobelActive™ 3.0	•						Unigrip	
NobelActive ASC Zirconia						•	Omnigrip	
Straumann						•	Star/Torx	
Zimmer®						•	0.050" (1.25 mm) hex	

Atlantis™ titanium and zirconia abutments utilize the same torque and driver setting as the original implant manufacturer.

Loaner Torque Kits Available Upon Request



Call: 519.742.2966 or 1.800.387.6533
 Email: info@silvadentallab.ca
 Hours: Monday to Friday
 8:30AM — 5:00PM